


Cell Injury Controller II



The **Cell Injury Controller II** is an electronically controlled pneumatic device that provides all the control and measurement functions necessary for producing strain-induced trauma to in vitro cell cultures. The CIC II facilitates the study of morphologic, physiologic and biochemical responses of cultured cells to strain-induced trauma.

Tissue cultures are grown in a commercially available tissue culture system from Flexcell Int'l. Corp. Controlled gas pressure is used to induce a physiologically rapid deformation of the tray's Silastic bottom resulting in quantifiable biaxial stretch of the cells adhered to the Silastic surface. Injury severity is determined by controlling the flow of gas in and out of the sealed culture well. The peak culture well pressure is captured providing an accurate indicator of the resulting biaxial stretch.

The CIC II accepts both the Flex I® 29.45cm² culturing trays (used with the earlier CIC Model 94A) and the BioFlex® 57.75cm² culturing trays. Depending on the cell type, the degree of injury, and the culture conditions the injured cells may die or repair. Therefore, the system can be employed to study the responses to trauma, including cell injury, repair, death or pharmacologic intervention.

See publications dealing with biaxial strain induced injury of brain derived cells in culture at:

<http://www.people.vcu.edu/~eellis/CICPubs.html>

For information or orders, contact:

Guy M. Blankinship III
Virginia Commonwealth University
Medical Center
1101 East Marshall Street, Room B1-060
P.O. Box 980470
Richmond, Virginia 23298-0470
Phone: (804) 828-1902
Fax: (804) 828-1903
E-mail: gblankin@vcu.edu



Flex I and BioFlex are registered trademarks of Flexcell International Corp., Hillsborough, N.C.

Silastic is a registered trademark of Down Corning Corp., Midland, Mich.

Teflon is a registered trademark of DuPont Corp., Wilmington, Del.